



Innovations in Molecular Separation

MAB Purification with SMB Chromatography and Semba Octave™ System



Semba Octave™ Chromatography System



- A versatile bench top 8-column system
- Capable of performing SMBC and other continuous automated separation protocols
- Suitable for milligram-to-multigram scale purification
- 100% biocompatible flow path
- SembaPro™ software controls multiple protocols
- Compatible with many commercially available purification columns and chromatography process media



Height: 27" / 69 cm


Width: 21.5" / 57 cm

Depth: 21" / 54 cm

Weight: 150 lb / 68 kg

MAb Purification - Protein A

- Parameters for Protein A SMBC
 - POROS[®] MabCapture[™] A (8 x 1.0 ml columns)
 - Step Mode Configuration
 - Feed flow rate 5.0 ml/min concentrated tissue culture fluid (1 g/L antibody concentration)
- Results
 - Productivity 200 mg/h
 - Purity >99%
 - Recovery 95%



MAb Purification - Single Column vs Step SMBC

- MAb were purified using single column affinity purification with a similar yield (95%) and purity > 95%

| Purification Mode | Protein A Resin | Buffer | Time |
|-------------------|-----------------|--------|-------|
| Single Column | 100 ml | 960 ml | 6.2 h |
| Step-SMBC | 8 x 1 ml | 320 ml | 1.8 h |

SMBC for Protein Purification

- Semba Octave™ System accelerates process development and scale-up times leading to gram amounts of pure protein faster
 - Two fractions to analyze: extract (target protein) and raffinate (untagged proteins, unwanted cellular components)
 - No product variability – fractions collected after the system reaches equilibrium do not change
 - Ideally suited for affinity purifications, SEC and other separations
 - The CIP step can be incorporated in the protocol controlled by the SembaPro™ software